


From Brain Signals to Cognition

EEG Workshop: Learning EEG Research Methods from Scratch

Time: 9 am – 1 pm, May 8, 2021 (Saturday) Mode: via ZOOM

By Dr Guang Ouyang, Assistant Professor, Academic Unit of Human Communication, Development, and Information Sciences, Faculty of Education, HKU



Cognition

Outline



- Basic concepts
- Major principles and methodologies
- Real-time demonstration of EEG data collection
- Data visualization and analysis
- Advanced topics


Feedback



A dose of history and philosophy



Hans Berger (1873 – 1941) discovered brain EEG and its association with cognition



EEG



Cognition?




Brain in a vat



What 100+ year of EEG research brings to us?



EEG



- Functional activities
- Cognitive activities
- Cognitive abilities
- Emotional states
- Mental states
- Pathological states
- Development/aging
- etc


Some recent fancy stuff

Brain-activity-to-text decoder




Makin, J. G., Moses, D. A., & Chang, E. F. (2020). Machine translation of cortical activity to text with an encoder–decoder framework. *Nature neuroscience*, 23(4), 575-582.

Levels of neural system



What generates EEG?



Electrical signals


Related to
cognitive activity

Electrical cable



What generates EEG?


- **Mechanism:** parallelly oriented neurons generate current activity, causing fluctuations of electrical signals on the scalp (EEG) and cortex (ECoG). The current further induces magnetic field (MEG).



Hari, R., & Parkkonen, L. (2015). The brain timewise: how timing shapes and supports brain function. *Philosophical Transactions of the Royal Society B: Biological Sciences*, 370(1668), 20140170.

An important note

- EEG/MEG directly measures neural activities.
- It possesses rich information in neural temporal dynamics (at milliseconds)
- which is not accessible by many other technologies (e.g., MRI, fNIRS) that measure hemodynamics



Hari, R., & Parkkonen, L. (2015). The brain timewise: how timing shapes and supports brain function. *Philosophical Transactions of the Royal Society B: Biological Sciences*, 370(1668), 20140170.

Basic features




Outline

- Basic concepts
- Major principles and methodologies
- Real-time demonstration of EEG data collection
- Data visualization and analysis
- Advanced topics



How to do cognitive research using EEG?



How to do cognitive research using EEG?

1. ERP method (Event-related potential)

Scenario 1:

stimulus

Evoked response
activity




hit

Scenario 2:

2. Spectrum analysis method



change of mental states

ERP method



ERP (Event-related Potential)


ERP method




ERP (Event-related Potential)

ERP method

Change in timing
(reflects cognitive processing speed)



Change in amplitude
(reflects cognitive effort)




ERP method

How ERP components reflect cognitive activity?

Early components

- Low level processing, e.g., size, luminance, sound intensity, etc
- Sometimes by top-down modulation, e.g., emotion, attention, etc.



Late components


- Medium to high level processing, e.g., conflict processing, surprises, emotion, memory recollection, etc.

ERP (Event-related Potential)


ERP method

Examples of stimulus manipulation

Luminance: **high**, **low**




Occurrence probability
of stimulus:
low, **high**




Tone/note differentiation
do do do do re do do

People with
good ability



People with
bad ability




Condition 1

Condition 2

Semantic violation

I like to eat coffee




ERP method

Why do I need to check ERP when I can already see relevant information in behavior?



1. You don't know yourself that well (subjective feeling can be imprecise)
2. Much information is not available in behavioral data

Test the your mental speed!




<https://humanbenchmark.com/tests/reactiontime>

Reaction Time



Spectrum analysis method


- Change of brain's internal state in a long-lasting way
- Examples:
 - Feeling excited/alerted/anxious/calm/drowsy...
 - Mood/emotion change
 - Cognitive load change
 - Meditation exercise
 - Different vegetative level
 - ...



2. Spectrum analysis method


Spectrum analysis method

Oscillation




Spectrum analysis method

Signal 1



10 Hz


Signal 2




30 Hz

1 second

Spectrum analysis method



Spectrum analysis method



Spectrum analysis method

Delta wave (1-3Hz)

Functional association

- Sleep wave
- mixed

Theta wave (4-7Hz)

Functional association

- Active control
- Attention

Alpha wave (8-12Hz)

Functional association

- Idle state
- Cognitive load

Beta wave (13-30Hz)

Functional association


- Voluntary movements

Gamma wave (>30Hz)

Function

Functional association

- Cognitive load
- Intense neural computation
- Increased sensorimotor activity




Outline

- Basic concepts
- Major principles and methodologies
- Real-time demonstration of EEG data collection
- Data visualization and analysis
- Advanced topics






Set-up of an EEG acquisition system



Some technical points worth mentioning

- Ambient noise, environmental signal
- Gel-based (wet) electrode, dry electrode, semi-dry electrode
- Lab EEG, portable EEG, Indoor & outdoor
- Some latest fancy technologies:
 - Unobtrusive EEG
 - Concealed EEG
 - Earphone/earbud EEG
 - Headset EEG
 - Connecting with smartphones





Product: **mBraintrain smarting system**
<https://mbraintrain.com/>


Outline

- Basic concepts
- Major principles and methodologies
- Real-time demonstration of EEG data collection
- Data visualization and analysis
- Advanced topics



1. ERP method (event-related potential)

Analyzing the brain's response to a specific event




Condition 1

Condition 2

2. Spectrum analysis method


Analyzing the change in the brain state during a long period of time



More convenient
Less flexible

Processing and analysis

UI based



Script based

```
clear;
eeglab nogui

%% data loading
data_path = 'D:\EEG_workshop\data\';
%need to change to your actual data path

part1 = 'sample_visual_oddball';
EEG = pop_loadset('filename',[part1,'.set'],'filepath',data_path);

%% pre-processing
%bandpass filtering
EEG = pop_eegfiltnew(EEG, 'locutoff',1,'hicutoff',45,'plotfreqs',0);
%detect bad channel and interpolate
std_temp = std(EEG.data,1,2);
ol = find(isoutlier(std_temp,'ThresholdFactor',4)&std_temp>mean(std_temp));
if ~isempty(ol) EEG = pop_interp(EEG, ol, 'spherical');end
%average referencing
EEG = pop_reref( EEG, []);
%run ica (the stop criterion usually needs to be lower)
EEG = pop_runica(EEG, 'isotype', 'runica', 'extended',1,'stop',0.001,'interrupt','on');
%automatically detect artifacts and clean the data
[comps,info] = MARA(EEG);
EEG = pop_subcomp( EEG, comps, 0);

%% generate ERP
epoch_twd = [-200,1000];%epoch time window
bl_twd = [-200, 0]; %baseline time window

ERPs = get_erp_func(EEG,epoch_twd,bl_twd,{'S 21'},{'S 22'});

ch = find(strcmpi([EEG.chanlocs.labels],'Pz'));
figure;plot(ERPs.t_axis,squeeze(ERPs.ERPs(ch,1,:)));
xlabel('time after stimulus (ms)');
ylabel('amplitude (\u00b5V)');
title('brain response pattern');
```


Raw EEG data



Less convenient
More flexible

Results

ERP



Understanding the mechanism of coding even if you don't code

```
clear;
eeglab nogui

%% data loading
data_path = 'D:\EEG_workshop\data\';
%need to change to your actual data path

parti = 'sample_visual_oddball';
EEG = pop_loadset('filename',[parti,'.set'],'filepath',data_path);


%% pre-processing

%bandpass filtering
EEG = pop_eegfiltnew(EEG, 'locutoff',1,'hicutoff',45,'plotfreqz',0);
%detect bad channel and interpolate
std_temp= std(EEG.data,1,2);
ol = find(isoutlier(std_temp,'ThresholdFactor',4)&std_temp>mean(std_temp));
if ~isempty(ol) EEG = pop_interp(EEG, ol, 'spherical');end
%average referencing
EEG = pop_reref( EEG, []);
%run ica (the stop criterion usually needs to be lower)
EEG = pop_runica(EEG, 'icatype', 'runica', 'extended',1,'stop',0.001,'interrupt','on');
%automatically detect artifacts and clean the data
[comps,info] = MARA(EEG);
EEG = pop_subcomp( EEG, comps, 0);

%% generate ERP
epoch_twd = [-200,1000];%epoch time window
bl_twd = [-200, 0]; %baseline time window

ERPs = get_ERP_func(EEG,epoch_twd,bl_twd,{'S 21'},{'S 22'});

ch = find(strcmpi({EEG.chanlocs.labels},'Pz'));
figure;plot(ERPs.t_axis,squeeze(ERPs.ERPs(ch,:,:)));
xlabel('time after stimulus (ms)');
ylabel('amplitude (\muV)');
title('brain response pattern');
```



General routine

Data visualization and analysis

To be covered in another
focused workshop

Raw data

Preprocessing

May include:

- Re-sampling
- Re-referencing
- Band-pass filtering
- Removing artifactual components
- etc





Visualization and Analysis

- ERP
- Spectrum
- Statistical analysis



Examples of artifacts (fake brain signals)


Typical:



Disconnected:





Drifting



Reason why we need to deal with artifacts

Condition 1 (task 1): speaking

Condition 2 (task 2): hearing a speech



Difference can be merely from the fact that the first one is speaking and the second one is not.



EEG analysis softwares

- Commercial softwares from EEG equipment companies
- Matlab + EEGLAB
- Matlab + Letswave
- Matlab + Fieldtrip
- Matlab + Brainstorm
- MNE-Python
- Many others


- 
1. ERP analysis
 2. Spectrum analysis

1. ERP analysis

Based on a visual oddball task:



average
→





ERP (Event-related Potential)





Code and materials:




Data visualization and analysis




Difference map:




100-150




178-228




256-306




333-383




411-461




489-539




567-617




644-694



722-772



800-850



Statistics

on large data sample


Individual events
Anecdotal findings
Special cases
Personal experiences
Intuitive feeling





Generalizable effects
Universal principle
Theories
Scientific evidences

Data visualization and analysis


From an individual subject:




Averaged from 30 subjects:



ERP amplitude difference (300-600ms) for all subjects



Prob. Dist.



Advanced statistical analyses


The EEG data collection and preprocessing will generate data that are not different from any other kinds of data in other fields

e.g., humanity and social sciences, biological sciences, psychology, educational sciences.


So, after processing and parameterizing your data, you can apply any kinds of statistical analysis depending on your experimental design.

- Correlation
- ANOVA/ANCOVA, etc
- Linear mixed model
- Factor analysis
- Structural equation modeling
- etc



2. Spectrum analysis







Sample data



Data visualization and analysis



Difference map:



Outline

- Basic concepts
- Major principles and methodologies
- Real-time demonstration of EEG data collection
- Data visualization and analysis
- Advanced topics



Upon acquiring the basic skill sets and understanding of EEG data and technology, you can further ...

- **Conduct sophisticated analyses**
 - Time-frequency analysis, network/connectivity analysis, extract complex features (e.g., entropy), dynamical modeling, advanced statistic modelling, etc
- **Study complex questions**
 - Real-life scenes, social interaction, teaching and learning, classroom setting, complex problem solving, complex decision making, etc

Give your feedback

